

AOX assay — well by well

Same two layouts used for both temperatures: fill A15 and A27 from the Plate A list, B15 and B27 from the Plate B list. Work down each column of the list in order. Three additions per well, always in this order.

PLATE A

well	1 — base	2 — azide slot, 50 µL	3 — inhibitor slot, 20 µL	final in vial
A1	5.03 mL cells	50 µL water	20 µL ethanol	—
A2	5.03 mL cells	50 µL azide 500 mM	20 µL SHAM 50 mM	azide 5 + SHAM 0.2
A3	5.03 mL cells	50 µL azide 500 mM	20 µL nPG 25 mM	azide 5 + nPG 0.1
A4	5.03 mL cells	50 µL azide 500 mM	20 µL ethanol	azide 5
A5	5.03 mL cells	50 µL azide 50 mM	20 µL ethanol	azide 0.5
A6	5.03 mL cells	50 µL azide 500 mM	20 µL SHAM 50 mM	azide 5 + SHAM 0.2
B1	5.03 mL cells	50 µL azide 500 mM	20 µL nPG 25 mM	azide 5 + nPG 0.1
B2	5.03 mL cells	50 µL water	20 µL ethanol	—
B3	5.03 mL cells	50 µL azide 200 mM	20 µL ethanol	azide 2
B4	5.03 mL cells	50 µL azide 500 mM	20 µL SHAM 50 mM	azide 5 + SHAM 0.2
B5	5.03 mL cells	50 µL azide 1 M neat	20 µL ethanol	azide 10
B6	5.03 mL cells	50 µL azide 500 mM	20 µL nPG 25 mM	azide 5 + nPG 0.1
C1	5.03 mL cells	50 µL azide 500 mM	20 µL ethanol	azide 5
C2	5.03 mL cells	50 µL azide 100 mM	20 µL ethanol	azide 1
C3	5.03 mL cells	50 µL azide 500 mM	20 µL nPG 25 mM	azide 5 + nPG 0.1
C4	5.03 mL PLAIN YMS (no cells)	50 µL water	20 µL ethanol	nothing
C5	5.03 mL cells	50 µL azide 500 mM	20 µL SHAM 50 mM	azide 5 + SHAM 0.2
C6	5.03 mL cells	50 µL water	20 µL ethanol	—
D1	5.03 mL cells	50 µL azide 50 mM	20 µL ethanol	azide 0.5
D2	5.03 mL cells	50 µL azide 1 M neat	20 µL ethanol	azide 10
D3	5.03 mL cells	50 µL azide 100 mM	20 µL ethanol	azide 1
D4	5.03 mL cells	50 µL azide 500 mM	20 µL ethanol	azide 5
D5	5.03 mL cells	50 µL azide 200 mM	20 µL ethanol	azide 2
D6	5.03 mL PLAIN YMS (no cells)	50 µL water	20 µL ethanol	nothing

PLATE B

well	1 — base	2 — azide slot, 50 µL	3 — inhibitor slot, 20 µL	final in vial
A1	5.03 mL cells	50 µL azide 200 mM	20 µL SHAM 50 mM	azide 2 + SHAM 0.2
A2	5.03 mL cells	50 µL azide 500 mM	20 µL nPG 250 mM	azide 5 + nPG 1.0
A3	5.03 mL cells	50 µL water	20 µL SHAM 50 mM	SHAM 0.2
A4	5.03 mL cells	50 µL azide 200 mM	20 µL nPG 25 mM	azide 2 + nPG 0.1
A5	5.03 mL cells	50 µL azide 500 mM	20 µL SHAM 250 mM	azide 5 + SHAM 1.0
A6	5.03 mL cells	50 µL water	20 µL nPG 250 mM	nPG 1.0
B1	5.03 mL cells	50 µL azide 500 mM	20 µL SHAM 250 mM	azide 5 + SHAM 1.0
B2	5.03 mL cells	50 µL azide 200 mM	20 µL SHAM 50 mM	azide 2 + SHAM 0.2
B3	5.03 mL cells	50 µL water	20 µL nPG 25 mM	nPG 0.1
B4	5.03 mL cells	50 µL azide 500 mM	20 µL ethanol	azide 5
B5	5.03 mL cells	50 µL azide 200 mM	20 µL nPG 25 mM	azide 2 + nPG 0.1
B6	5.03 mL PLAIN YMS (no cells)	50 µL water	20 µL ethanol	nothing
C1	5.03 mL cells	50 µL azide 200 mM	20 µL nPG 25 mM	azide 2 + nPG 0.1
C2	5.03 mL cells	50 µL water	20 µL SHAM 250 mM	SHAM 1.0
C3	5.03 mL cells	50 µL azide 500 mM	20 µL nPG 250 mM	azide 5 + nPG 1.0
C4	5.03 mL cells	50 µL azide 200 mM	20 µL SHAM 50 mM	azide 2 + SHAM 0.2
C5	5.03 mL cells	50 µL water	20 µL nPG 250 mM	nPG 1.0
C6	5.03 mL cells	50 µL azide 500 mM	20 µL SHAM 250 mM	azide 5 + SHAM 1.0
D1	5.03 mL cells	50 µL azide 500 mM	20 µL ethanol	azide 5
D2	5.03 mL cells	50 µL water	20 µL SHAM 50 mM	SHAM 0.2
D3	5.03 mL cells	50 µL azide 500 mM	20 µL nPG 250 mM	azide 5 + nPG 1.0

D4	5.03 mL cells	50 µL water	20 µL nPG 25 mM	nPG 0.1
D5	5.03 mL PLAIN YMS (no cells)	50 µL water	20 µL ethanol	nothing
D6	5.03 mL cells	50 µL water	20 µL SHAM 250 mM	SHAM 1.0

Every vial is 5.03 mL + 50 µL + 20 µL = 5.10 mL. Ethanol lands at 0.39% and water at 0.98% in all 24, including the blanks. Fill all four plates alternately from the one 500 mL cell batch, swirling between draws. The blanks are the only vials without cells. Azide: never acidify, never down the sink.